

Medline Rubber Drain Tube (8.0 x 2.0mm)

Product Code: 211494

Brand: Medline Scientific

The Medline Rubber Drain Tube (8.0 x 2.0mm) is manufactured from 45° natural rubber, high strength, good resistance to tear, fatigue, and abrasion. Has excellent electrical properties but poor resistance to ozone, mineral, and petroleum based oils.



Specifications:

Temperature	-65°C to +70°C
Applications	Transferring gases, liquids, and diluted acids
Coil Length (m)	20

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com