

## Medline Silicone Tube (3.0 x 0.75mm)


**Product Code:** 21300415

**Brand:** Medline Scientific

With excellent resistance to temperature, ozone and weathering, the Medline Silicone Tube (3.0 x 0.75mm) offers increased tear strength, less bacterial build up, and a higher clarity over peroxide silicone rubber tubing.

It also has no smell, taste taint or reaction with chlorine and is less likely to have surface extrusion lines, which are sometimes visible on peroxide silicone tubing.

Manufactured from FDA approved materials.



### Specifications:

<b>Temperature</b>	-40°C to +200°C (+250°C short term)
<b>Applications</b>	Food, medical, and pharmaceutical. Peristaltic pump use.
<b>Coil Length (m)</b>	15

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)