

Medline Silicone Tube (8.0 x 2.0mm)


Product Code: 21301715

Brand: Medline Scientific

With excellent resistance to temperature, ozone and weathering, the Medline Silicone Tube (8.0 x 2.0mm) offers increased tear strength, less bacterial build up, and a higher clarity over peroxide silicone rubber tubing.

It also has no smell, taste taint or reaction with chlorine and is less likely to have surface extrusion lines, which are sometimes visible on peroxide silicone tubing.

Manufactured from FDA approved materials.



Specifications:

Temperature	-40°C to +200°C (+250°C short term)
Applications	Food, medical, and pharmaceutical. Peristaltic pump use.
Coil Length (m)	15

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com