

## Medline Six Place Multi-Position Extraction Mantle (250ml)

**Product Code:** EAM-9202-06

**Brand:** Medline Scientific

Designed for distillation and extraction, Medline's Six Place Multi-Position Extraction Mantle (250ml) includes support rods and scaffolding.

The body is made of steel with a chemical resistant coating and the cover of high grade stainless steel.



### Specifications:

<b>Volume (ml)</b>	250
<b>Positions</b>	6

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)