

Medline Small PET Sweet Jars With Caps (1500ml)

Product Code: 870004

Brand: Medline Scientific

Ideal for the safe and convenient disposal of bench top waste and other small contaminated waste in microbiology or pathology, Medline's PET Sweet Jars With Caps are made of high clarity polyethylene terephthalate making them ideal to autoclave.

Specifications:

Volume (ml)	1500
Height (mm)	150
Diameter (mm)	125
Neck Diameter (mm)	120
Material	Polyethylene Terephthalate
Sterility	Non-Sterile
Quantity	30



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com