

Medline Small PET Sweet Jars With Caps (1500ml)

Product Code: 870004

Brand: Medline Scientific

Ideal for the safe and convenient disposal of bench top waste and other small contaminated waste in microbiology or pathology, Medline's Small PET Sweet Jars With Caps (1500ml) are made of high clarity polyethylene terephthalate.

Specifications:

Volume (ml)	1500
Height (mm)	150
Diameter (mm)	125
Neck Diameter (mm)	120
Material	Polyethylene Terephthalate
Sterility	Non-Sterile
Quantity	30

