

## Medline Ultra-Flat Magnetic Stirrer with Auto Reverse Function

**Product Code:** FEI02001D

**Brand:** Medline Scientific

The Medline Ultra-Flat Magnetic Stirrer with Auto Reverse Function is very thin and lightweight and can be set to mix and stir up to 1000ml between 15 and 1500rpm. To deliver an excellent mixing performance, this stirrer can carry out automatic reverse rotations every 30 seconds. Its distinctive set-up plate and casing is chemically resistant and durable.

The low-profile and smooth design of the Medline Ultra-Flat Magnetic Stirrer allows it to be used in compact spaces such as fume hoods or cabinets. It is ideal in teaching laboratories as it can be stacked in cupboards when not in use to save space. The top-plate is virtually indestructible, waterproof and very easy to clean. It can also be used in cold rooms and incubators (non CO<sub>2</sub>).



### Features

- OLED display provides easy reading of the speed and time
- Touch key is more user-friendly and durable
- Special black tempered glass set-up plate is waterproof, shatterproof; very easy to clean and dissipate heat
- The most advanced magnetic coil technology allows for wear-free drive without moving parts
- Automatic reverse rotation every 30 seconds provides excellent mixing performance
- Max. stirring quantity (H<sub>2</sub>O) of 1000ml
- Speed range of 15 - 1500rpm

## Specifications

<b>Working Plate Material</b>	Tempered Glass
<b>Speed Range (rpm)</b>	15 - 1500
<b>Mixing Action</b>	Selectable Auto-Reverse
<b>Display</b>	OLED
<b>Max. Stirring Quantity (H<sub>2</sub>O)</b>	1000ml
<b>Max. Stirring Bar (mm)</b>	25
<b>Stirring Surface (mm)</b>	Ø 100
<b>Input Power</b>	6W
<b>Dimensions (W x D x H)</b>	117 x 12 x 180mm
<b>Weight (kg)</b>	0.3

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)