

Medline Urine Collection Kit (Boric Acid) - Screw Cap Primary Tube with Boric Acid & Collection Cup with Collection Cup

Product Code: 40906BKITR

Brand: Medline Scientific

The Medline Urine Collection Kit (Boric Acid) - Screw Cap Primary Tube with Boric Acid & Collection Cup is designed for convenience of use in hospitals, GP surgeries, and health clinics.

Boric acid acts as a preservative and maintains microbiological quality of samples while also preventing overgrowth of microorganisms.

Boxes are divided into 9 bags of 50 kits, making distribution to GPs or wards easy. Kits consist of primary tube and collection cup with spout. Screw cap compatible with major analysers.

