

## Medline Vortex Mixer (VM-3000MD)

**Product Code:** VM-3000MD

**Brand:** Medline Scientific

Equipped with a strong and silent motor, Medline's Vortex Mixer (VM-3000MD) produces a controlled and reproducible agitation throughout the speed range utilizing either touch activation or continuous modes.

Suction cups displayed on its base are designed to minimise vibration and movement.

### Specifications:

<b>Speed (rpm)</b>	20 to 3000
<b>Dimensions (mm)</b>	130 x 155 x 155
<b>Weight (kg)</b>	3.1
<b>Power</b>	220V - 50/60Hz



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)