

Medline Vortex Mixer (VX-3000MD)

Product Code: VX-3000MD

Brand: Medline Scientific

The Medline Vortex Mixer (VX-3000MD) features a strong, reliable and silent motor that produces regulated and efficient agitation throughout the speed range utilising either touch activation or continuous mode. An LED display is included for optimum control and display of time and speed. Rubber feet on the base minimise any vibration and movement.



Specifications

Shaking Motion	Orbital
Orbital Diameter (mm)	6
Speed Range (rpm)	0 to 3000
Operating Mode	Touch, Continuous
Dimensions; W x D x H (mm)	107 x 146 x 166
Weight (kg)	3

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com