

## Hettich 220 - 240V Refrigerated Benchtop Centrifuge (MIKRO 200R)

**Product Code:** HETGB-C2205

**Brand:** Hettich

The Hettich 220 - 240V Refrigerated Benchtop Centrifuge (MIKRO 200R) are suitable for research configuration tasks (e.g. PCR, DNA & RNA). Its Fast Cool function refrigerates quickly (down to + 4°C in 10 to 15 minutes), making the MIKRO 200 perfect for any fast paced laboratory. Quick and efficient refrigeration guarantees that even temperature sensitive samples are centrifuged correctly.



### Features

- Robust and steady, temperature control system with pre-cool function from -10 to + 40°C (optional)
- For DNA, RNA, and temperature sensitive samples
- Bio-contamination, autoclavable, phenol resistant lids
- **Fast Cool Function:** Down to + 4°C in 10 - 15 minutes
- Run-Silent™ Operation (≤ 58dB, softer than normal conversation)
- Easy-Lift™/Safety-Lock™ lid
- Five year warranty
- FDA registered centrifuge manufacturer

### Specifications

<b>Max. Capacity</b>	30 x 1.5/2.0ml
<b>Max. RPM/RCF</b>	15,000/21,382
<b>Temperature Control</b>	-10 to + 40°C
<b>Dimensions (H x W x D)</b>	10.23" x 11" x 13.5"
<b>Weight</b>	61.72lbs

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

Tel: (+44) 01865 400321

Email: [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

Website: [www.medlinescientific.com](http://www.medlinescientific.com)