

Medline 100 to 1500rpm Analogue Stirrer (MS-200)

Product Code: MS-200

Brand: Medline Scientific

Constructed from chemically resistant ceramic coated steel, Medline's 100 to 1500rpm Analogue Stirrer (MS-200) delivers excellent reliability with reproducible results at a sensible price. A standard 12.5mm screwed thread is included for the addition of a retort rod (not included).



Specifications

Description	Stirrer
Plate Area (mm)	180 x 180
Speed Range (rpm)	100 to 1500
Dimensions (W x D x H)	410 x 210 x 150mm
Weight (kg)	4

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com