

Medline 100 to 1500rpm Analogue Hotplate Stirrer (MS-300)

Product Code: MS-300

Brand: Medline Scientific

Medline's 100 to 1500rpm Analogue Hotplate Stirrer (MS-300) offers a large heating area of 180 x 180mm with safe controlled heating up to 400°C. Constructed from chemically resistant ceramic coated steel, this lab hotplate stirrer delivers excellent reliability with reproducible results at a sensible price. A screwed 12.5mm screwed thread is included for the addition of a retort rod (not included).



Specifications

Description	Hotplate Stirrer
Plate Area (mm)	180 x 180
Temperature Range (°C)	Max. 400
Speed Range (rpm)	100 to 1500
Dimensions (w x d x h)	410 x 210 x 150mm
Weight (kg)	5

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com