


## Deltalab Stainless Steel Rack for 24 Containers (58mm)

**Product Code:** N-240

**Brand:** Deltalab

Maximise space for your sample containers with Deltalab's Stainless Steel Rack for 24 Containers (58mm). This wire constructed stainless steel rack can hold up to 24 containers of 58mm Ø and has 63mm square holes.



### Specifications:

<b>Dimensions (l x w x h)</b>	393 x 262 x 45mm
<b>Containers Capacity</b>	24 (6 x 4)
<b>Colour</b>	White
<b>Material</b>	Stainless Steel
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)