

Jeio Tech 150 Litre Economy Forced Convection Oven (OF-21E)

Product Code: I4032

Brand: Jeio Tech

The Jeio Tech 150 Litre Forced Convection Oven (OF-21E) is an affordable, multi-functional fast drying laboratory oven with vertical airflow. It is corrosion resistant and has a convenient presetting function for the three most commonly temperature settings.

Other features include microprocessor PID control, adjustable dual slide vents for inner vapour circulation, and easy-to-grip circular door handles.

Features:

Performance:

- Temperature range from ambient +10°C to 220°C.
- Microprocessor PID control/temperature calibration/automatic tuning.
- Forced air convection generated by a circular fan provides high and uniform drying performance.
- Adjustable dual slide vents for inner vapour circulation.

Convenience:

- Intuitive control panel with bright LED displays (0.1°C resolution).
- Dual wait on/off timer modes (1 min. to 99 hr. 59 min.).
- Convenient presetting function for three most commonly used temperature settings.
- Corrosion-resistant incoloy sheath heating element and stainless steel interior.
- Easy-to-grip circular door handle.

Safety:

- Keypad lock function prevents accidental changes during operation.
- Automatic run after power interruption.
- Over-temperature and over-current protection.
- Audible and visible open door alarms.



Specifications:

Chamber Volume (litres)	150	
Temperature Range (°C)	Amb. +10 to 220	
Temperature Fluctuation at 100°C	± 0.4	
Temperature Variation at 100°C	± 2.3	
Dimensions	Internal (w x d x h)	585 x 460 x 555mm
	External (w x d x h)	740 x 625 x 922mm
	No. of Wire Shelves (standard/max.)	2/14
	Net Weight (kg)	59
Electrical Requirements (230V, 50/60Hz)	7.4A	

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com