

RBS 50 Detergent (5 Litres)

Product Code: RBS50

Brand: RBS

RBS 50 Detergent (5 Litres) is a non-foaming liquid alkaline detergent used in laboratories for cleaning glassware and laboratory equipment made of stainless steel, plastic, porcelain, quartz, and ceramic.

It also removes organic and mineral residues, proteins, chemical solutions, distillation residues, fats, oils, and cell cultures.

In the industrial sector, it is used as a technical product for high performance cleaning of metal parts, optical lenses, glass, and electronic components.

RBS 50 complies with EN 1040, EN 1276, EN 1275, and NFT 72-190 (spectrum-4) standards under certain conditions and concentrations.

Contains chlorine.

Specifications

Size (litres)	5 (bottle)
Use Concentration	1% v/v (normal/manual cleaning) and 0.3 to 0.5% v/v (ultrasonic bath)



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com