

## RBS A 155 Cleaning Agent (5 Litres)

**Product Code:** RBSA155

**Brand:** RBS

RBS A 155 Cleaning Agent (5 Litres) is a non-foaming cleaning and decontaminating agent used for glassware and laboratory equipment in washing machines. This is thanks to its high content in active alkaline elements and non-foaming property.

It can also be used for surfaces made of stainless steel, glass, plastic, and allows total removal of soils (e.g. fats, oils, chemical solutions, blood and distillation residues).

RBS A 155 is chlorine and phosphates-free.



## Specifications:

<b>Size (litres)</b>	5 (bottle)
<b>Use Concentration</b>	0.3 to 2% v/v (normal cleaning)

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)