

## RBS A 375 Neutraliser & Acid Rinse (5 Litres)

**Product Code:** RBSA375

**Brand:** RBS

RBS A 375 Neutraliser & Acid Rinse (5 Litres) is suitable for pre-wash and rinsing cycles in automatic washing machines.

It is used in the first rinsing cycle after the alkaline washing with either RBS A 155 (liquid), RBS 50 pF (liquid), or RBS SOLID pF (powder).

RBS A 375 allows the neutralisation and the removal of residual traces of carbonates, hydroxides, and calcareous elements on the cleaned items and on the internal components of the washing machine.

It can also be used during the pre-wash cycle when glassware and materials to be cleaned are soiled with minerals.

It also does not contain surfactants or phosphorus that may interfere with certain types of analysis (e.g. water and phosphates) or pollutes the environment.



### Specifications:

<b>Size (litres)</b>	5 (bottle)
<b>Use Concentration</b>	0.2 to 0.5% v/v depending on the type of machine and the soils to be removed

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)