

RBST 230 Detergent (5 Litres)

Product Code: RBST230

Brand: RBS

RBST 230 Detergent (5 Litres) is a multi-purpose neutral cleaning agent for sensitive laboratory equipment, soaking, and ultrasonic baths.

Thanks to its high content in active ingredients and its neutral pH, RBST 230 is the ideal product for cleaning sensitive glassware and materials to bases or acids. It presents excellent compatibility with most of the materials.

At the recommended working concentration, the RBST 230 is neither corrosive, irritating, nor toxic for both the cleaned items and the user.



Specifications

Size (litres)	5 (bottle)
Use Concentration	4 to 6% v/v in water

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com