

RBST 305 Detergent (5 Litres)

Product Code: RBST305

Brand: RBS

RBST 305 Detergent (5 Litres) is a strong acidic detergent used for manual cleaning, degreasing, and de-scaling surfaces and laboratory equipment.

It is also recommended for cleaning animal cages and rooms in laboratories and veterinary institutions.



Specifications

Size (litres)	5 (bottle)
Use Concentration	20 to 50ml of RBST 305 per litre of water

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com