


Deltalab Red Caps for Transparent Polypropylene Screw Thread Microtubes

Product Code: 409112/4

Brand: Deltalab

Made of medical grade polypropylene, Deltalab's Red Caps for Transparent Polypropylene Screw Thread Microtubes feature an internal O-ring to ensure a leakproof seal.



Specifications

Colour	Red
Dimensions (mm)	13 x 6
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com