

Hettich 208 - 240V Refrigerated Benchtop Centrifuge, Without Rotor (ROTINA 420R)

Product Code: MD5211156

Brand: Hettich

The Hettich 208 - 240V Refrigerated Benchtop Centrifuge, Without Rotor (ROTINA 420R) is a high-performance benchtop centrifuge for large volumes, yet not overly large in size. This new, compact benchtop centrifuge is for high sample throughput as well as for high sample volumes. With a maximum capacity of 4 x 600ml, the ROTINA 420R is suitable for applications in clinical chemistry, biotechnology, and life science research facilities.

Features

- Robust and steady, temperature control system with pre-cool function from -20 to + 40°C
- Quiet operation
- Easy-Lift™/Safety-Lock™ lid
- Five year warranty

Specifications

Max. Capacity	4 x 600ml
Max. RPM/RCF	15,000/24,400
Temperature Control	-20 to + 40°C
Dimensions (H x W x D)	16.6" x 28" x 25.7"
Weight	257.9lbs



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com