

## Jeio Tech 35 Litre Cubic Non-Vacuum Desiccator (SDC-30)

**Product Code:** AAAD5011

**Brand:** Jeio Tech

The stackable Jeio Tech 35 Litre Cubic Non-Vacuum Desiccator (SDC-30) is made up of a moulded cubic design, providing space efficiency. This desiccator is shatterproof, non-corrosive, and has an easy-to-clean polycarbonate construction.

The internal space is maximised by using a maximum of three shelves and its space saving, stackable design allows efficient use of cramped lab space. Humidity levels are checked using the hygrometer attached to the front door.



### Features

- One moulded design enables virtually no leakage and exceptional durability
- High-quality silicone gasket ensures a greaseless, airtight seal for a clean and cost effective test environment
- Highly transparent, shatterproof, non-corrosive, and easy to clean polycarbonate construction
- Maximised utilisation of internal space using a maximum of three shelves
- Stackable space saving design allowing efficient use of cramped lab space
- Easy humidity level checking with the hygrometer attached to the front door

## Specifications

<b>Internal Volume (L)</b>	35
<b>Max. Load per Shelf or Drawer (kg)</b>	5
<b>Number of Shelves (included/max.)</b>	1/3
<b>Number of Drawers (included/max.)</b>	1/3
<b>Internal Dimensions (W x D x H)</b>	322 x 367 x 302mm
<b>External Dimensions (W x D x H)</b>	360 x 410 x 326mm
<b>Weight (kg)</b>	4.9

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)