

Jeio Tech 35 Litre Cubic Non-Vacuum Desiccator (SDC-30U)

Product Code: AAAD5111

Brand: Jeio Tech

The stackable Jeio Tech 35 Litre Cubic Non-Vacuum Desiccator (SDC-30U) is made up of a moulded cubic design, providing space efficiency. This desiccator is shatterproof, non-corrosive, and has an easy-to-clean polycarbonate construction.

The internal space is maximised by using a maximum of three shelves and its space saving, stackable design allows efficient use of cramped lab space. Humidity levels are checked using the hygrometer attached to the front door. This laboratory desiccator also features a UV blocking polycarbonate material to minimise the damage to, or discolouration of light-sensitive samples.



Features

- One moulded design enables virtually no leakage and exceptional durability
- High-quality silicone gasket ensures a greaseless, airtight seal for a clean and cost effective test environment
- Highly transparent, shatterproof, non-corrosive, and easy to clean polycarbonate construction
- Maximised utilisation of internal space using a maximum of three shelves
- Stackable space saving design allowing efficient use of cramped lab space
- Easy humidity level checking with the hygrometer attached to the front door
- Minimising damages or discolouration of light-sensitive samples

Specifications

Internal Volume (L)	35
Max. Load per Shelf or Drawer (kg)	5
Number of Shelves (included/max.)	1/3
Number of Drawers (included/max.)	1/3
Internal Dimensions (W x D x H)	322 x 367 x 302mm
External Dimensions (W x D x H)	360 x 410 x 326mm
Weight (kg)	4.9

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com