

Jeio Tech 50 Litre Cubic Non-Vacuum Desiccator (SDC-45)

Product Code: AAAD5021

Brand: Jeio Tech

The stackable Jeio Tech 50 Litre Cubic Non-Vacuum Desiccator (SDC-45) is made up of a moulded cubic design, providing space efficiency. This desiccator is shatterproof, non-corrosive, and has an easy-to-clean polycarbonate construction.

The internal space is maximised by using a maximum of five shelves and its space saving, stackable design allows efficient use of cramped lab space. Humidity levels are checked using the hygrometer attached to the front door.



Features

- One moulded design enables virtually no leakage and exceptional durability
- High-quality silicone gasket ensures a greaseless, airtight seal for a clean and cost effective test environment
- Highly transparent, shatterproof, non-corrosive, and easy to clean polycarbonate construction
- Maximised utilisation of internal space using a maximum of five shelves
- Stackable space saving design allowing efficient use of cramped lab space
- Easy humidity level checking with the hygrometer attached to the front door

Specifications

Internal Volume (L)	50
Max. Load per Shelf or Drawer (kg)	5
Number of Shelves (included/max.)	2/5
Number of Drawers (included/max.)	1/5
Internal Dimensions (W x D x H)	322 x 367 x 422mm
External Dimensions (W x D x H)	360 x 410 x 447mm
Weight (kg)	6.2

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com