

Ceti Triple Scale Hand-Held Digital Refractometer (SR-95)

Product Code: 8950.0000M

Brand: Ceti

The CETI Triple Scale Hand-Held Digital Refractometer (SR-95) has a triple range of Brix, Refractive Index, and Salinity coupled with the convenience of an RS232 interface. Other conditions include automatic temperature compensation, Centigrade and Fahrenheit measurements, a stainless steel sample plate, and instantaneous digital results.



Specifications

Type(s)	Brix	Refractive Index	Salinity
Range	0 - 95%	1.3330 - 1.5318%	0 - 30%
Resolution	0.1%	0.0001%	0.1%
Accuracy	+/- 0.1%	+/- 0.00015%	+/- 0.1%

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com