

CETI Stereo-I0P Student Microscope

Product Code: I111.0010M

Brand: CETI

The CETI Stereo-I0P Student Microscope is a high-quality, robust microscope ideal for educational purposes and allows students to comfortably view specimens with its binocular vertical optical head.

This built-to-last stereo microscope has a super-wide 20mm field of view with built-in dioptric compensation on one eye piece, providing simple parfocality adjustment. It comes standard with 2x Plan Achromatic objectives.



Specifications:

Optical Head	Binocular, vertical	
	Interpupillary Distance Adjustment (Sidentopf Type): 51 to 75mm	
	Diotropic compensation on one eyepiece	
Eye Pieces	10x/20mm wide-field	
Objectives	2x Planachromatic (supplied as standard)	
Stage	60mm Ø Frosted and black/white contrast plate	
	Two stage clips	
Focusing Device	Coarse focus clamp of optical head onto retort stand	
	Fine focus	
	A pair of knobs on each side of the stand	
	Rack and pinion system	
Packing	Dimensions (w x d x h)	250 x 115 x 170mm
	Gross Weight (kg)	1.3
	Additional Packaging Contents	Eyecups and dust cover

Unit 3, Tower Business Park
 Warpsgrove Lane
 Chalgrove
 Oxfordshire
 OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com