

## CETI Stereo-30 Student Microscope

**Product Code:** I I I I.0030M

**Brand:** CETI

Specially designed for education purposes, the robust CETI Stereo-30 Student Microscope can be used for various laboratory applications and delivers an outstanding optical performance, allowing students to carry out long productive working sessions.

This stereo microscope is equipped with a binocular optical head (inclined at 45° and 360° - rotatable) and can be operated with 2x or 4x Plan Achromatic objectives.



### Specifications:

<b>Optical Head</b>	Binocular, inclined at 45° and 360° (rotatable)
	Interpupillary Distance Adjustment: 55 to 75mm
<b>Eyepieces</b>	10x/20mm wide field
<b>Objectives</b>	2x or 4x Plan Achromatic
<b>Stage</b>	95mm Ø Frosted and black/white contrast plate
	Two stage clips
<b>Illumination</b>	Variable light intensity
	Episcopic 12V/10W Tungsten
	Diascopic 12V/10W Halogen
<b>Mains Power</b>	220 - 240V/50 - 60Hz

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)