

Trajan Series I 90° Frosted Plain Edged Microscope Slides (Double Frosted)

Product Code: TJ471022100

Brand: Trajan

Trajan's Series I 90° Frosted Plain Edged Microscope Slides (Double Frosted) are non-coated and frosted at one end, making it easier for laboratories to identify samples. These slides are designed for printing with automated slide printers.

All slides are pre-cleaned and inspected before distribution.

Features

- Manufactured using high-quality German extra plane float glass
- Packed with double cellophane in box-type cardboard boxes

Specifications

Colour	None (Double Frosted)
Dimensions (mm)	76 x 26
Thickness (mm)	1
Quantity	50 per pack

