

Trajan Series I 90° Frosted Ground Edged Microscope Slides (Double Frosted)

Product Code: TJ471042100

Brand: Trajan

Trajan's Series I 90° Frosted Ground Edged Microscope Slides (Double Frosted) are non-coated and frosted at one end, making it easier for laboratories to identify samples. These slides are designed for printing with automated slide printers.

All slides are pre-cleaned and inspected before distribution.

Features

- Manufactured using high-quality German extra plane float glass
- Packed with double cellophane in small cardboard boxes

Specifications

Colour	None (Double Frosted)
Dimensions (mm)	76 x 26
Thickness (mm)	1
Quantity	50 per pack



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com