

Trajan Pink Series I 90° Frosted Ground Edged Microscope Slides

Product Code: TJ471042203

Brand: Trajan

Trajan's Pink Series I 90° Frosted Ground Edged Microscope Slides are non-coated and frosted at one end, making it easier for laboratories to identify samples. These slides are designed for printing with automated slide printers.

All slides are pre-cleaned and inspected before distribution.

Features

- Manufactured using high-quality German extra plane float glass
- Packed with double cellophane in small cardboard boxes

Specifications

| | |
|------------------------|-------------|
| Colour | Pink |
| Dimensions (mm) | 76 x 26 |
| Thickness (mm) | 1 |
| Quantity | 50 per pack |



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com