

## Trajan White Series I 90° Adhesive Ground Edged Microscope Slides

**Product Code:** TJ47104222I

**Brand:** Trajan

Trajan's White Series I 90° Adhesive Ground Edged Microscope Slides have a specially designed coating that produces a good charged surface for histology, cytology, and microbiology applications providing an excellent and consistent performance. The clear charge promotes cell and tissue adhesion, and provide a better base for aqueous reagents.

These slides are suitable for frozen sectioning, heat-induced antigen/epitope retrieval (HIER) immunohistochemistry (IHC) and in situ hybridization (ISH) based techniques.

All slides are pre-cleaned and inspected before distribution.

### Features

- Manufactured using high-quality German extra plane float glass
- Packed with double cellophane in lid-type cardboard boxes

### Specifications

<b>Colour</b>	White
<b>Dimensions (mm)</b>	76 x 26
<b>Thickness (mm)</b>	1
<b>Quantity</b>	72 per pack



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)