

Trajan White Series | 90° Advanced Adhesive Ground Edged Microscope Slides

Product Code: TJ471042241

Brand: Trajan

Trajan's White Series | 90° Advanced Adhesive Ground Edged Microscope Slides have a hydrophilic surface with a contact angle around 30° and strong adherence.

These slides are suitable for diagnostic applications including immunohistochemistry (IHC) and can withstand hard operations. With auto fluorescence, they are also ideal for in situ hybridization (ISH) based applications.

All slides are pre-cleaned and inspected before distribution.

Features

- Manufactured using high-quality German extra plane float glass
- Packed with double cellophane in lid-type cardboard boxes

Specifications

Colour	White
Dimensions (mm)	76 x 26
Thickness (mm)	1
Quantity	72 per pack



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com