

## Trajan White Series 2 Adhesive Large Frost Microscope Slides

**Product Code:** TJ472042491

**Brand:** Trajan

Trajan's White Series 2 Adhesive Large Frost Microscope Slides are manufactured using fine glass with excellent flatness and transmittance with minimal intrinsic fluorescence.


All slides are carefully cut, pre-cleaned, and inspected to ensure that they meet your laboratory's needs.

### Features

- Optimized with immunostaining with minimal background staining and adhesion
- Slides are optimized for use with a range of automated immunostainers
- Manufactured using high quality extra plane float glass
- Packed with double cellophane in lid-type cardboard boxes

### Specifications

<b>Colour</b>	White
<b>Frost Size (mm)</b>	20
<b>Dimensions (mm)</b>	76 x 26
<b>Quantity</b>	100 per pack



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)