

Trajan Yellow Series 3 Adhesive Large Frost Microscope Slides

Product Code: TJ473042492

Brand: Trajan

Trajan's Yellow Series 3 Adhesive Large Frost Microscope Slides are manufactured with a specialized fortified adhesion coating that provides improved tissue sample adherence and extensibility. Each slide is made of fine glass with excellent flatness and transmittance with minimal intrinsic fluorescence.

All slides are carefully cut, pre-cleaned, and inspected to ensure that they meet your laboratory's needs.

Features

- Optimized for immunostaining with fortified adhesion
- Slides are optimized for use with a range of automated immunostainers
- Manufactured using high quality extra plane float glass
- Packed with double cellophane in lid-type cardboard boxes

Specifications

Colour	Yellow
Frost Size (mm)	20
Dimensions (mm)	76 x 26
Quantity	100 per pack



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com