

Trajan Series I No. 1 Circular Coverslips (20mm)

Product Code: TJ471112000

Brand: Trajan

Trajan's Series I No. 1 Circular Coverslips (20mm) are manufactured using ultra-thin pure white borosilicate glass in controlled conditions with uniform thickness and flatness in non-corrosive and non-fogging packaging.

Their special non-stick treatment prevents the coverslips sticking together. Each coverslip can be chosen together with Trajan's microscope slides for both manual and automated coverslipping.



Specifications:

Diameter (mm)	20
Material	Glass
Sterility	Non-Sterile
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com