

## Trajan Series I No. 1 Circular Coverslips (5mm)

**Product Code:** TJ471110600

**Brand:** Trajan

Trajan's Series I No. 1 Circular Coverslips (5mm) are manufactured using ultra-thin pure white borosilicate glass in controlled conditions with uniform thickness and flatness in non-corrosive and non-fogging packaging.

Their special non-stick treatment prevents the coverslips sticking together. Each coverslip can be chosen together with Trajan's microscope slides for both manual and automated coverslipping.



### Specifications:

<b>Diameter (mm)</b>	5
<b>Material</b>	Glass
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)