

Jeio Tech 2 Litre Ultrasonic Cleaner (UCP-02)

Product Code: BDH532011K

Brand: Jeio Tech

Operating at a cleaning frequency of 40kHz, the Jeio Tech 2 Litre Ultrasonic Cleaner (UCP-02) provides maximum cleaning power and strong reliability.

Perfect for cleaning components in different sizes, its #304 stainless steel tank is easy-to-clean and sealed with a durable/flame resistant cover. This cleaner also features a simple digital LED display and temperature controller.



Features:

- Durable/flame resistance ABS material used and cover with #304 stainless steel tanks.
- High-frequency 40 kHz sound waves provide greater cleaning power and increased reliability.

Specifications:

Bath Volume (litres)		2
Temperature Range (°C)		-
Ultrasonic Power (W)		100
Dimensions	Bath Interior (w x l, d)	142 x 132, 104mm
	Exterior (w x l x h)	290 x 210 x 200mm
	Net Weight (kg)	2
Electrical Requirements (230V, 50/60Hz)		0.6A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com