

Jeio Tech 10 Litre Ultrasonic Cleaner (UCP-10)

Product Code: BDH532012K

Brand: Jeio Tech

Operating at a cleaning frequency of 40kHz, the Jeio Tech 10 Litre Ultrasonic Cleaner (UCP-10) provides maximum cleaning power and strong reliability.

Perfect for cleaning components in different sizes, its #304 stainless steel tank is easy-to-clean and sealed with a durable/flame resistant cover. This cleaner also features a simple digital LED display and temperature controller.

Features:

- Durable/flame resistance ABS material used and cover with #304 stainless steel tanks.
- Over-temperature protection with alarm.
- Digital LED display using splash-proof keypads: timer (1 to 99 min., max. 4 hour operating) and temperature monitor.
- High-frequency 40 kHz sound waves provide greater cleaning power and increased reliability.
- Multiple output mode (Normal/Pulse/Sweep) with power selection (Low/Middle/High).
- Before cleaning, dissolved gas removed by degas function.
- Drain valve, basket, cover as standard offer.



Specifications:

Bath Volume (litres)		10
Temperature Range (°C)		Amb. to 70
Ultrasonic Power (W)		300
Dimensions	Bath Interior (w x l, d)	295 x 240, 160mm
	Exterior (w x l x h)	490 x 368 x 297mm
	Net Weight (kg)	11
Electrical Requirements (230V, 50/60Hz)		1.3A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com