

## Jeio Tech 20 Litre Ultrasonic Cleaner (UCP-20)

**Product Code:** BDH532013K

**Brand:** Jeio Tech

Operating at a cleaning frequency of 40kHz, the Jeio Tech 20 Litre Ultrasonic Cleaner (UCP-20) provides maximum cleaning power and strong reliability.

Perfect for cleaning components in different sizes, its #304 stainless steel tank is easy-to-clean and sealed with a durable/flame resistant cover. This cleaner also features a simple digital LED display and temperature controller.

### Features:

- Durable/flame resistance ABS material used and cover with #304 stainless steel tanks.
- Over-temperature protection with alarm.
- Digital LED display using splash-proof keypads: timer (1 to 99 min., max. 4 hour operating) and temperature monitor.
- High-frequency 40 kHz sound waves provide greater cleaning power and increased reliability.
- Multiple output mode (Normal/Pulse/Sweep) with power selection (Low/Middle/High).
- Before cleaning, dissolved gas removed by degas function.
- Drain valve, basket, cover as standard offer.



## Specifications:

<b>Bath Volume (litres)</b>		20
<b>Temperature Range (°C)</b>		Amb. to 70
<b>Ultrasonic Power (W)</b>		500
<b>Dimensions</b>	<b>Bath Interior (w x l, d)</b>	500 x 290, 160mm
	<b>Exterior (w x l x h)</b>	490 x 368 x 297mm
	<b>Net Weight (kg)</b>	17
<b>Electrical Requirements (230V, 50/60Hz)</b>		2.2A

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)