

Jeio Tech 10 Litre Ultrasonic Cleaner (UCS-10)

Product Code: BDH532032K

Brand: Jeio Tech

Providing state-of-the-art ultrasonic cleaning, the Jeio Tech 10 Litre Ultrasonic Cleaner (UCS-10) is user-friendly with a clear digital LED display, containing splash-proof keypads.

Its tank is made of durable #304 stainless steel and has various output modes.

Features:

- External and tank parts are made of durable 304 stainless steel.
- Over-temperature protection with alarm.
- Digital LED display using splash-proof keypads: timer (1 to 99 min., max. 4 hours operating), temperature monitor.
- High-frequency 40 kHz sound waves provide greater cleaning power and increased reliability.
- Multiple output mode (normal/pulse/sweep) with power selection (low/middle/high).
- Drain valve, basket, cover as standard offer.



Specifications:

Bath Volume (litres)		10
Temperature Range (°C)		Amb. to 70
Ultrasonic Power (W)		300
Dimensions	Bath Interior (w x l, d)	295 x 240, 160mm
	Exterior (w x l x h)	325 x 330 x 320mm
	Net Weight (kg)	11
Electrical Requirements (230V, 50/60Hz)		1.3A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com