

Arctiko ULTF 283 Litre -86°C Ultra Low Temperature Chest Freezer

Product Code: ULTF 320

Brand: Arctiko

The Arctiko ULTF 283 Litre -86°C Ultra Low Temperature Chest Freezer shows that high performance freezing with focus on low energy consumption is the future. This ultra low temperature freezer features Arctiko's advanced microprocessor controller, which ensures that all control, alarm and data logging features are included as standard. In addition, it is equipped with Arctiko's original single compressor technology.

Features

- Low energy consumption with insulated inner doors
- Heated door frame for frost-free gasket and easy door opening
- 100% HCFC/CFC free refrigerants and low noise level
- Microprocessor controller with digital display of temperature, battery level, ambient temperature and temperature graph
- Battery backup (approx. 72 hours) for alarms, data, logging and temperature display in case of power outage



Specifications

Dimensions	Internal (W x D x H)	1100 x 440 x 630mm
	External (W x D x H)	1262 x 698 x 885mm
Capacity (L)		283
Noise level (dB)		<50
Temperature range (°C)		-86/-40
Max. ambient temperature (°C)		20
Power supply (V)		230
Frequency (Hz)		50/60
Power consumption (kWh/24h)		7.8

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com