

## Arctiko ULUF 7 Litre -86°C Ultra Low Temperature Freezer

**Product Code:** ULUF 15

**Brand:** Arctiko

Arctiko's ULUF 7 Litre -86°C Ultra Low Temperature Freezer will add a new dimension to your laboratory. This laboratory freezer can be used as a bench top or under bench model. Its compact space-saving design and low noise level makes it an excellent storage solution for your laboratory. It is produced with Arctiko's original single compressor technology.

### Features:

- Low energy consumption with insulated inner doors
- Heated door frame for frost-free gasket and easy door opening
- 100% HCFC/CFC free refrigerants and low noise level
- Microprocessor controller with digital display of temperature, battery level, ambient temperature and temperature graph
- Battery backup (approx. 72 hours) for alarms, data, logging and temperature display in case of power outage.



## Specifications

<b>Dimensions</b>	<b>Internal (W x D x H)</b>	150 x 143 x 310mm
	<b>External (W x D x H)</b>	400 x 630 x 665mm
<b>Capacity (L)</b>		7
<b>Noise level (dB)</b>		<55
<b>Temperature range (°C)</b>		-86/-40
<b>Max. ambient temperature (°C)</b>		25
<b>Power supply (V)</b>		230
<b>Frequency (Hz)</b>		50/60
<b>Power consumption (kWh/24h)</b>		4.6

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)