

Arctiko ULUF 642 Litre -86°C Ultra Low Temperature Freezer

Product Code: ULUF 750-2M

Brand: Arctiko

Arctiko's ULUF 642 Litre -86°C Ultra Low Temperature Freezer will add a new dimension to your laboratory. This laboratory freezer can be used as a bench top or under bench model. Its compact space-saving design and low noise level makes it an excellent storage solution for your laboratory. It is produced with Arctiko's original single compressor technology.

Features:

- Low energy consumption with insulated inner doors
- Heated door frame for frost-free gasket and easy door opening
- 100% HCFC/CFC free refrigerants and low noise level
- Microprocessor controller with digital display of temperature, battery level, ambient temperature and temperature graph
- Battery backup (approx. 72 hours) for alarms, data, logging and temperature display in case of power outage.



Specifications

Dimensions	Internal (W x D x H)	790 x 680 x 1345mm
	External (W x D x H)	1030 x 885 x 1990mm
Capacity (L)		642
Noise level (dB)		< 55
Temperature range (°C)		-86/-40
Max. ambient temperature (°C)		25
Power supply (V)		230
Frequency (Hz)		50/60
Power consumption (kWh/24h)		16.9

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com