

## Arctiko 610 Litre -90°C Integra Ultra Low Temperature Freezer

**Product Code:** ULUF P610

**Brand:** Arctiko

The Arctiko 610 Litre -90°C Integra Ultra Low Temperature Freezer has a dual cooling system, which consists of two individual systems that present excellent security. If one of the cooling systems fails, the other system will still control the temperature at -80°C. Therefore, your samples are unlikely to be affected by a rise in temperature. This feature, provided for by Arctiko's original and highly reliable single compressor technology, could be the difference between losing or saving a freezer full of priceless samples.

### Features

- Dual, independent cooling system for double security
- Maximum cooling of compressors with airflow control and filter-less technology
- VIP panels and insulated inner doors to ensure lowest energy consumption and ECO mode to further reduce energy consumption
- Reduced noise level down to 49dB.
- Microprocessor control with clear display and self-diagnostic functions
- Adjustable shelves and easily removable inner doors
- 100% HCFC/CFC free



## Specifications

Dimensions	Internal (W x D x H)	746 x 600 x 1320mm
	External (W x D x H)	1035 x 860 x 1999mm
Capacity (L)		610
Noise level (dB)		<49
Temperature range (°C)		-90/-40
Max. ambient temperature (°C)		25
Power supply (V)		230
Frequency (Hz)		50/60
Power consumption (kWh/24h)		17.1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)