

Jeio Tech 16 Litre Cylindrical Vacuum Desiccator (VDP-25G)

Product Code: AAAD3111

Brand: Jeio Tech

With a high-quality silicone O-ring and specially designed 'Locking Clip', Jeio Tech's 16 Litre Cylindrical Vacuum Desiccator (VDP-25G) offers light weight yet flexible strength and provides outstanding vacuum capability with virtually no leakage.

Features:

- Maximum vacuum of 1 Torr (1.33×10^{-4} Mpa) for more than 72 hours.
- High-quality silicone gasket ensures greaseless airtight seal for a clean test environment.
- Highly transparent, shatterproof, non-corrosive, and easy to clean polycarbonate construction.
- Three-way stopcock provides great convenience for consistent and uniform vacuum draw, vacuum release or gas exchange.
- A flange lock holds the top and the bottom together during non-vacuum conditions.
- A vacuum gauge is attachable at the centre of the lid for easy vacuum level checking.
- The top and bottom of the desiccator are compatible with each other for diverse use.
- Ideal for storing tall equipment, gas substitution, general drying, and preserving.



Specifications:

Internal Volume (litres)		16
Gauge Range (Mpa)		0 to 0.1
Nozzle (mm)		9.5
Max. Permissible Vacuum		1.33×10^{-4} Mpa (Analogue)
Material	Lid/Bottom	Polycarbonate
	Body/Base	
	Perforated Sample Tray	Polypropylene
	Drying Agent Tray/Locker	
	Vacuum Seal	Silicone Rubber
Dimensions	Exterior (i.d. x h)	308 x 374mm
	Net Weight (kg)	2.7

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com