

Jeio Tech 6 Litre Round Vacuum Desiccator (VDR-20G)

Product Code: AAAD2021

Brand: Jeio Tech

With a high-quality silicone O-ring and specially designed 'Locking Clip', the Jeio Tech 6 Litre Round Vacuum Desiccator (VDR-20G) offers light weight yet high tensile strength and provides outstanding leakproof vacuum capability.

The built-in vacuum gauge on this laboratory desiccator allows fast and accessible analysis of the vacuum level.

Features:

- Maximum vacuum of 1 Torr (1.33×10^{-4} Mpa) for more than 72 hours.
- High-quality silicone gasket ensures greaseless air-tight seal for a clean test environment.
- Highly transparent, shatterproof, non-corrosive, and easy-to-clean polycarbonate construction.
- Three-way stopcock provides great convenience for consistent and uniform vacuum draw, vacuum release, or gas exchange. The barbed fitting on each side secures attachment of tubing (maximum 2 three-way stopcocks are installable if required when placing an order).
- Flange lock holding the top and bottom together during non-vacuum conditions.



Specifications:

Internal Volume (litres)		6
Gauge Range (Mpa)		0 to -0.1
Nozzle (mm)		9.5
Max. Permissible Vacuum (bodies)		1.33×10^{-4} Mpa
Material(s)	Lid/Bottom	Polycarbonate
	Body/Base for VDP	
	Perforated Sample Tray	Polypropylene
	Drying Agent Tray	
	Locker	
	Vacuum Seal	Silicone Rubber
Dimensions	Exterior (i.d. x h)	Ø 242 x 279mm
	Net Weight (kg)	1.5

Unit 3, Tower Business Park
 Warpsgrove Lane
 Chalgrove
 Oxfordshire
 OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com