

## Jeio Tech 20 Litre Round Vacuum Desiccator (VDR-30G)

**Product Code:** AAAD2221

**Brand:** Jeio Tech

With a high-quality silicone O-ring and specially designed 'Locking Clip', the Jeio Tech 20 Litre Round Vacuum Desiccator (VDR-30G) offers light weight yet high tensile strength and provides outstanding leakproof vacuum capability.

The built-in vacuum gauge on this laboratory desiccator allows fast and accessible analysis of the vacuum level.



### Features:

- Maximum vacuum of 1 Torr ( $1.33 \times 10^{-4}$  Mpa) for more than 72 hours.
- High-quality silicone gasket ensures greaseless air-tight seal for a clean test environment.
- Highly transparent, shatterproof, non-corrosive, and easy-to-clean polycarbonate construction.
- Three-way stopcock provides great convenience for consistent and uniform vacuum draw, vacuum release, or gas exchange. The barbed fitting on each side secures attachment of tubing (maximum 2 three-way stopcocks are installable if required when placing an order).
- Flange lock holding the top and bottom together during non-vacuum conditions.

## Specifications:

<b>Internal Volume (litres)</b>		20
<b>Gauge Range (Mpa)</b>		0 to -0.1
<b>Nozzle (mm)</b>		9.5
<b>Max. Permissible Vacuum (bodies)</b>		$1.33 \times 10^{-4}$ Mpa
<b>Material(s)</b>	<b>Lid/Bottom</b>	Polycarbonate
	<b>Body/Base for VDP</b>	
	<b>Perforated Sample Tray</b>	Polypropylene
	<b>Drying Agent Tray</b>	
	<b>Locker</b>	
	<b>Vacuum Seal</b>	Silicone Rubber
<b>Dimensions</b>	<b>Exterior (i.d. x h)</b>	Ø 385 x 399mm
	<b>Net Weight (kg)</b>	4.2

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)