

## Deltalab Square Holed Test Tube Rack with Grip (for 13mm Tubes)

**Product Code:** W-011

**Brand:** Deltalab

Autoclavable and stackable when unloaded, Deltalab's Square Holed Test Tube Rack with Grip (for 13mm Tubes) can be dismantled for space saving.

Its unique grip feature ensures that the tubes are firmly held in the rack. An alphanumeric grid reference allows you to easily locate samples. This rack is very easy to assemble.

### Specifications:

<b>Max. Diameter of Tubes (mm)</b>	13
<b>Tube Capacity</b>	6 x 15 (90)
<b>Colour</b>	White
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)