

TECHNICAL DATA SHEET

Standard flocked swab BP 80 mm

code 310206K

code	description	presentation	breakpoint mm	case quantity	volume case
310206K	standard flocked swab	flow-pack	80 mm	2 x 1000	0.029

Directive 93/42/CEE. Medical Devices.



INTENDED USE- DESCRIPTION

Sterile invasive system for collecting samples by direct contact with the patient in a surgical environment. Indicated for the collection and transport of biological samples for subsequent microbiological analysis.

Standard flocked swab with polystyrene backing and flocked head, with excellent absorption and elution of microorganisms.

Characteristics :



- Support made of polystyrene --> Optimum safety.
- Head swab made of flocked fibers--> increased absorption and elution of the sample on the media.
- Maximum adaptability, wide range.
- Maximum comfort: anatomical and ergonomic design.
- Breakpoint at 80 mm.
- Flow-pack packaging.

Swab: sterile by radiation (R).

Designed for cell culture and compatible with molecular diagnostic techniques.

Swab dimensions

Length: 152 mm.



WRITTEN	VERIFIED	AUTHORISED
Joan Ibáñez Marketing Department	Anna Mir Technical, Development, Quality and Environment Director	Mónica Torras Sales and Marketing Director